

ACCUVIN

D- & L-Lactic Acid for Tomato




Cat. No. 615 or 620

Intended Use

The D- & L-Lactic Acid Kit is intended for measuring the combined D-Lactic Acid plus L-Lactic Acid level of tomato products as an indicator of Lactic Acid Bacteria contamination.

Methodology

D- & L-Lactic Acid is based on the change in color exhibited by a tetrazolium color indicator during a reaction involving Lactic acid and nicotine-adenine dinucleotide (NAD) in the presence of the enzymes D-Lactate dehydrogenase and L-Lactate dehydrogenase.



Sample

Samples of tomato juice can be used as they are. Samples of tomato paste, ketchup/catsup and sauce should be diluted 1:5 prior to analysis using the Quick Dilute 1:5 tubes. The ACCUVIN D- & L-Lactic Acid patent pending test strip removes the usual interferences from colored and turbid samples. Samples do not have to be pre-filtered or treated with color removing substances such as activated carbon or polyamide powder. Sample temperature may be from 5°C - 35°C (41°F - 95°F).

Procedure

1. Squeeze upper sampler bulb. Dip sampler tip into tomato juice or diluted tomato product sample, then release sampler bulb to aspirate sample. (If you prefer to use an air displacement pipette, set sample volume at 55 µL.)
2. Transfer sample to the rectangular absorbent layer on back of test strip by squeezing sampler bulb. **Apply slight pressure with sampler tip for about 5 seconds while dispensing sample.** Note that only sample present in the sampler tip will be dispensed. Sample in the center, overflow bulb will not be dispensed. Wait 4 min. for color development.
3. Determine sample D- & L-Lactic Acid level in mg/L by comparing the developed color to the color chart on the test strip container. If test strip color falls between two color chips read numbers on second row. Note that if a sample was diluted prior to analysis, the sample lactic acid level is 5 times the level obtained from the color chart. **(Since fluorescent lights have a green cast, color matching is best under incandescent or natural lighting.)**

Storage

Store away from direct sunlight at temperatures below 80°F. Keep dry. Product is satisfactory until the date printed on the test strip package label.

www.ACCUVIN.com

Seller's limitation of liabilities: Every effort has been made to ensure the material contained in this informational insert and the results obtained with AV test strips are as accurate as possible, but no warranty or fitness is implied. Buyer shall not in any event be entitled to, and seller shall not be liable for, indirect, special, incidental or consequential damages of any nature including, without being limited to, loss of profit, promotional or manufacturing expenses, overhead, injury to reputation or loss of customers. Buyers recovery from seller for any claim shall not exceed buyer's purchase price for the products irrespective of the nature of the claim, whether in contract, tort, warranty, or otherwise.

615B, June 2007

Summary Interpretation

The tomato is a unique fruit that provides an essential source of vitamin C, potassium and antioxidants, primarily lycopene. Approximately one third of the world's annual production of about 80,000,000 metric tons of tomatoes is processed into products such as tomato paste, tomato pulp (diced tomatoes), and tomato sauce.¹

Lactic acid bacteria (LAB), especially of the genera *Lactobacillus* and *Leuconostoc*,^{2, 3, 6} are the bacterial agents most commonly involved in the spoilage of tomato products. They have a high tolerance to the low pH values (in the range of 3.5 to 4.6⁴) found in tomato products. When present at high levels or for long periods they can unacceptably affect a product's organoleptic properties^{5, 7, 8, 9, 10} through the production of lactic acid and acetic acid, diacetyl, ethanol or other chemical components.

There is a need for a simple, rapid, accurate method for measuring unacceptable contamination by LAB. The traditional method requires 72 hours, and does not provide a positive result if a previously-contaminated batch has been reprocessed.¹¹ Monitoring for changes in titratable acidity caused by the contaminating bacteria does not offer satisfactory sensitivity. Measuring for diacetyl and acetoin has been proposed,^{2, 3} but not all LAB produce these metabolites at acceptable levels, and some other LAB species first produce them but later metabolize them, so there is no remaining evidence of ongoing contamination.¹² Analysis for uracil has also been proposed,¹¹ but this method requires expensive HPLC instrumentation plus there is frequently insufficient distinction between levels present in non-contaminated products and levels produced by contaminating LAB. Since all LAB produce D- and/or L-Lactic acid, monitoring for L-Lactic Acid has been suggested,¹⁴ but measuring L-Lactic Acid plus D-Lactic Acid is the method of choice.^{5, 15, 17, 19, 20}

Uncontaminated tomato fruit and tomato products contain low levels of D + L-Lactic acid, in the 0 – 200 mg/kg range.^{5, 13, 14, 15} It is considered indicative of contamination when the level of these acids reaches or exceeds 300 mg/kg^{5, 15, 16, 17} or 400 mg/kg.¹⁸

References

1. USDA and USDA Foreign Agricultural Service reports.
2. C. Luster, "A rapid and sensitive sterility monitoring technique for aseptically processed bulk tomato paste," *J. Food Sci.*, 43: 1046 – 1052 **1978**.
3. B. J. Juven and H. Weisslowicz, "Chemical changes in tomato juices caused by lactic acid bacteria," *J. Food Sci.*, 46: 1543 – 1545 **1981**.
4. W. A. Hayes, P. G. Smith and A. E. J. Morris, "The production and quality of tomato concentrates," *Crit. Rev. Food Sci. Nutr.*, 38 (7): 537 – 564 **1998**.
5. E. Vicini, S. Gherardi, A. Bonici, C. Zoni and A. Trifiro, "Sviluppo di metaboliti e alterazioni chimiche causate da alcuni tipi di microrganismi in derivati vegetali," *Industria Conserve*, 63: 39 - 44 **1988**.
6. K. J. Bjorkroth and H. J. Korkella, "Lactobacillus fructivorans spoilage of tomato ketchup," *J. Food Protect.*, 60 (5): 505 -509 **1996**.
7. S. Porretta, G. Poli, M. Rondelli and G. Zanangeli, "Quality evaluation of tomato pulp," *Food Chem.*, 47: 379 - 386 **1993**.
8. S. Porretta, G. Poli, M. Rondelli and G. Zanangeli, "Definition of tomato pulp quality using a multivariate QDA method," *Lebensmittel Wissenschaft und Technologie*, 25 (5): 487 - 490 **1992**.
9. S. Porretta, "Composition of commercial tomato pulp and its relationship to some typical sensory characteristics," *J. Food Comp. Anal.*, 6 (2): 185 - 194 **1993**.
10. S. Porretta, "Analysis of sensory and physicochemical data on commercial tomato puree with pattern recognition techniques," *Z. Lebensmittel-Untersuchung u –Forschung*, 197 (6): 531 - 536 **1993**.
11. A. Hidalgo, C. Pompei, A. Galli and S. Cazzola, "Uracil as an index of Lactic acid bacteria contamination of tomato products," *J. Agric. Food Chem.*, 53: 349 - 353 **2005**.
12. C. G. Mallidis and F. Samaris, "A simple, rapid and sensitive method for monitoring contamination of aseptically processed tomato paste," *Int. J. Food Sci. Technol.*, 22: 59 - 66 **1987**.
13. S. S. El Miladi, W. A. Gould and R. L. Clements, "Heat processing effect on starch, sugars, proteins, amino acids, and organic acids of tomato juice," *Food Technol.*, 23 (5): 93 - 95 **1969**.
14. K. Kriz, I. Kraft, M. Krook and D. Kriz, "Amperometric determination of L-Lactate based on entrapment of lactate oxidase on a biosensor: application to tomato paste and baby food," *J. Agric. Food Chem.*, 50 (12): 3419 - 3424 **2002**.
15. S. Porretta and E. Vicini, "Changes in tomato pulp quality caused by lactic acid bacteria," *Int. J. Food Sci. Technol.*, 28: 611 - 616 **1993**.
16. S. Hanewinkel-Meshkini and W. Hackmann, "Beitrag zur beurteilung von tomatenmark," *Deutsche Lebensmittel Rundschau*, 85 (11): 351 -353 **1989**.
17. F. Mazzei, A. Azzoni, B. Cavalieri, F. Botre and C. Botre, "A multi-enzyme bioelectrode for the rapid determination of total lactate concentration in tomatoes, tomato juice and tomato paste." *Food Chem.*, 55 (4): 413 - 418 **1996**.
18. anon., "Disciplinare di produzione della Denominazione di Origine Proteta "Pomodoro San Marzano dell'Agro Sarnese-Nocerino," Reg. CE n. 1263/96.
19. M. R. Ackland, E. R. Trehwella, J. E. Reeder and P. G. Bean, "The detection of microbial spoilage in canned food using TLC," *J. Appl. Bact.* 51: 277 - 281 **1981**.
20. C. Consiglieri, "Valutazione del contenuto di acido D- e L-lattico ed acetico in concentrati di pomodoro cold-break," *Industrie Alimentari*, 30 (296): 738 - 740 **1991**.